

Organoid Processing Protocol

Optimized for airway organoids cultured as in Hynds et al. *Methods in Molecular Biology* 2016:

https://link.springer.com/protocol/10.1007%2F7651_2016_5

Materials/Equipment:

15 ml falcon tubes

V-bottomed 96-well FACS plate

1.5 ml Eppendorf tubes

4% PFA

70% EtOH

PBS

[HistoGel](#) (1 ml aliquots in 1.5 ml Eppendorf tubes)

Tray of ice

Plate-compatible centrifuge

Heat block

Procedure:

1. Pipette 6 ml ice-cold 4% PFA into 15 ml falcon tubes (1 per condition).
2. Transfer organoids to 15 ml falcon tubes on ice using a 200 ul pipette with the ends of the tips cut (combining 6-12 wells of a 96-well plate for histology processing).
3. Fix organoids for 30 mins on ice. During this time, pre-cool a labelled V-bottomed 96-well plate on ice, set centrifuge to 4°C and warm an aliquot of HistoGel to 65°C on a hot block.
4. Centrifuge 15 ml falcon tubes at 400 x g for 5 mins to pellet organoids (pellet should be visible).
5. Remove PFA being careful not to disturb pellet.
6. Resuspend pellet in 200 ul ice-cold PBS and transfer to a pre-cooled, pre-labelled 96-well FACS plate using a 200 ul pipette with the ends of the tips cut off (1 well per condition).
7. Centrifuge plate at 400 x g for 5 mins.
8. Carefully remove supernatant from wells using a 1 ml pipette.
9. Resuspend the pellet of organoids in 120 ul HistoGel using a 200 ul pipette with the ends of the tips cut off.
10. Allow HistoGel to set on ice for 10 mins.
11. Remove HistoGel from wells by gently pushing a 200 ul pipette tip down the wall of the well.
12. Store in labelled 1.5 ml Eppendorf tubes containing 500 ul 70% EtOH until processing for histology using routine dehydration to paraffin methods.